

Quarterly Pipette Calibration (gravimetric)

Date _____
Temp. (C) _____

Analyst _____
Z-Factor _____

Balance _____
Serial # _____

Pipette #	
VOLUME (ml)	WEIGHT (g)
MEAN	
CORR. MEAN	
STD DEV	
% CV	
% INACC.	
PASS / FAIL ?	

Pipette #	
VOLUME (ml)	WEIGHT (g)
MEAN	
CORR. MEAN	
STD DEV	
% CV	
% INACC.	
PASS / FAIL ?	

Pipette #	
VOLUME (ml)	WEIGHT (g)
MEAN	
CORR. MEAN	
STD DEV	
% CV	
% INACC.	
PASS / FAIL ?	

Pipette #	
VOLUME (ml)	WEIGHT (g)
MEAN	
CORR. MEAN	
STD DEV	
% CV	
% INACC.	
PASS / FAIL ?	

Pipette #	
VOLUME (ml)	WEIGHT (g)
MEAN	
CORR. MEAN	
STD DEV	
% CV	
% INACC.	
PASS / FAIL ?	

Pipette #	
VOLUME (ml)	WEIGHT (g)
MEAN	
CORR. MEAN	
STD DEV	
% CV	
% INACC.	
PASS / FAIL ?	

Corrective Action Due To Failures? _____

Acceptance Criteria

% Inaccuracy must be less than 2.00

No replicate may be greater than 2% from true volume.

%CV must be less than 1.00