

TSS Benchsheet [02NOV15]

For use with eight hour minimum oven dry time

Facility:	
Method Number:	SM 2540D-1997

Set Up		Take Out	
Analyst:		Analyst:	
Date in the Oven		Date out of the oven	
Time Samples in the Oven		Time samples out oven	
Oven Temperture In (°C)		Oven Temperature Out (°C)	

Control Limits			
Minimum sample weight if less than 500mL of sample if used (g)	0.0010	Minimum oven drying time	8 hrs
Maximum sample weight (g)	0.2000		

Sample ID	Filter or Pan Number	Sample Volume Used (mL)	Sample Plus Filter (g)	Filter Only (g)	Residue Mass (g)	Result not valid if error in either column	
						TSS (mg/L)	Reported Value (mg/L)
COMMENTS							